



DISINFECTING CLEANER

**HEALTHCARE HARD-SURFACE
DETERGENT-DISINFECTANT**

ACTIVE INGREDIENTS

10% QAC blend; Surfactant

PRODUCT CHARACTERISTICS

Super-concentrated; Unfragranced;
Non-corrosive; Non-toxic; Non-staining
Non-caustic; Biodegradable

RECOMMENDED USES

Manual or high pressure cleaning &
disinfection of hard surfaces.
Medical Device cleaning & disinfection.
Eradication of mould, mildew and fungal
growth within healthcare facilities.

PPE REQUIREMENTS

**Latex or nitrile gloves and chemical
splash goggles** are recommended when
handling or diluting the concentrated
product.



PRODUCT DESCRIPTION

INSTRUCTIONS FOR USE

Shake well before dilution. Dilute **SteriTech Disinfecting Cleaner** as per required efficacy from the Product Testing table (below right).

A 1% dilution is 10mL of product with 1L of water, (or 2 dosing pumps with 5L of water), a 2% dilution is 20mL of product with 1L of water, (or 4 dosing pumps with 5L of water).

Shake well before use. Spray, pour or mop the diluted solution onto the surface, and wipe down, or mop thoroughly. Allow the minimum contact time of your required efficacy, as per the Product Testing table (below right) before rinsing, or using the surface again.

Dilute only enough for immediate use.

Dispose as per the instructions provided on the label

Available sizes:

25L, 5L Dosing Pump, 1L, 1L Dosing Bottle, 500mL

Recommended Dilution Factor:

1:100 or 1% (10mL per L of water)

Physical Form:

Concentrate: Burgundy liquid

Once diluted: Translucent pink-red liquid

Odour: Mild chemical smell

pH (Concentrate): 11-13 **pH (Diluted):** 8.5-9.5

Risk Classification: B / IIa

GMDN Code: 47631

GMDN Descriptor:

Medical device disinfection agent.

NRCS Registration:

NRCS/8054/242613/224

SAHPRA Medical Device Registration:

Awaiting Class Registration Call-up

Incompatibilities:

Anionic detergents, Natural Soaps, Strong oxidisers



EFFICACY TESTING:



	Tested To:	Time: (mins)	Dilution:	Conditions:
Bactericidal Efficacy	SANS 51276 (EN 1276): Enterococcus hirae, Escherichia coli, Pseudomonas aeruginosa & Staphylococcus aureus SANS 53727 (EN 13727): Enterococcus hirae, Pseudomonas aeruginosa & Staphylococcus aureus	5	1%	Dirty
Fungicidal & Yeasticidal Efficacy	SANS 51650 (EN 1650): Aspergillus brasiliensis & Candida albicans SANS 53624 (EN 13624): Aspergillus brasiliensis & Candida albicans	15	2% 1%	Dirty
Sporocidal Efficacy	SANS 53704 (EN 13704): Bacillus subtilis	15	1%	Dirty
Corrosivity	SANS 636 (In-house): Non-corrosive up to 7 hours. Do not exceed a 7 hour exposure time.			
Dermal irritation	OECD 439: Category 2 irritant.			
Cytotoxicity	ISO 10993-5: Cytotoxic potential. Do not apply to broken skin, open wounds or cell-culture containers.			