



INSTRUMENT MILK

Lubricating instrument soak

PRODUCT DESCRIPTION

Available sizes: 5L Dosing Pump, 1L Dosing Bottle

Dilution Factor: 10% (100mL per L)

Physical Form: Opaque White liquid

Odour: Mild chemical odour

pH of Concentrate: 7-8

pH of diluted solution: 5.0-8.0

Incompatibility: Do not mix with other chemicals

PPE REQUIREMENTS

Latex or nitrile gloves and chemical splash goggles, however OSH also requires good ventilation when mixing, handling or discarding chemicals.



For additional information about this product, or to give product feedback to the manufacturer, please scan this QR code.



INSTRUCTIONS FOR USE

For manual soaking: Dilute with water at 1:10 or 10% [100mL, or 4 pumps with a dosing pump, of PrimeSafe Instrument Milk per 1L of water], immediately prior to use. Pour the solution into an appropriate HDPE, glass or stainless steel receptacle. Dip pre-cleaned, rinsed & dried instruments into the solution prior to packaging and or sterilisation. Diluted solution may be used for 7 days whereafter it should be discarded and a fresh solution made up.

For automated instrument washers: Dilute with water at 1:50 or 2% [2mL per 1L of water], and administer according to the washer manufacturer's recommended guidelines.