

# STAIN REMOVER

*Instrument Stain Remover & Descaler*

## PRODUCT DESCRIPTION

**Available sizes:** 5L, 1L

**Dilution Factor:** 12.5% (125mL per L)

**Physical Form:** Translucent clear liquid

**Odour:** Mild chemical odour

**Incompatibility:** Not for use on plated or non-stainless steel instruments

**pH:** 2 - 2.4 (Once diluted: 2.5-3.5).

## PPE REQUIREMENTS

**Latex or nitrile gloves and chemical splash goggles**, however OSH also requires the use of good ventilation when mixing, handling or discarding chemicals.



For additional information about this product, or to give product feedback to the manufacturer, please scan this QR code.



## INSTRUCTIONS FOR USE

**For manual use:** Dilute with hot water (>50°C) at 12.5% [125mL, or 5 pumps with a dosing pump, of PrimeSafe Stain Remover with 1L of water], immediately prior to use. Pour the solution into an appropriate HDPE, glass or Stainless Steel receptacle. Place pre-cleaned, rinsed & dried instruments into the solution for 15 minutes. Remove from the solution and rinse thoroughly under running water prior to packaging and or sterilisation.

Single use solution – discard immediately after use.

**For ultrasonic baths and or automated instrument washers:** Add 125mL, or 5 pumps with a dosing pump, of PrimeSafe Stain Remover per 1L of water directly into the detergent compartment of the machine. Run the machine as per the manufacturer's instructions. Upon removing the instruments from the machine, the instruments should be thoroughly rinsed under running water prior to proceeding to the following step in the reprocessing cycle of the institution.