



SteriTech
SPECIALIST INFECTION PROTECTION

CLINICAL CLEANER



HEALTHCARE HARD-SURFACE DETERGENT-DISINFECTANT

ACTIVE INGREDIENTS

5% Isopropylol; 1.5% QAC blend;
0.5% Triamine; 0.75% Surfactant

PRODUCT CHARACTERISTICS

Ready-to-Use; Unfragranced;
Non-corrosive; Non-toxic; Non-staining
Non-caustic; Biodegradable

RECOMMENDED USES

Manual cleaning & disinfection of
critical hard surfaces and glassware.
Medical Device cleaning & disinfection.
Eradication of mould, mildew and fungal
growth within healthcare facilities.

PPE REQUIREMENTS

**Latex or nitrile gloves and chemical
splash goggles (for liquid only)** are
recommended.



For additional information about
this product, or to give product
feedback to the manufacturer,
please scan this QR code.

PRODUCT DESCRIPTION

Available sizes: 25L, 5L, 100mL Mist Spray;
500mL Foaming Trigger Spray,
300s Wipes Bucket (Wipe Size: 170mm x 170mm)
100s Wipes Tub (Wipe Size: 170mm x 140mm)

Dilution Factor (Spray): Ready to use

Material Specification (Wipes):

20-30 gsm polypropylene / viscose blend

Physical Form (Spray): Translucent liquid

Physical Form (Wipe): White, nonwoven wipe

Odour: Mild chemical smell

pH: 9-10

Risk Classification: B / II

GMDN Code: 47631

GMDN Descriptor:

Medical Device Disinfectant

NRCS Reg. no.:

ACT5GNR529/242613/040/0658

Medical Device Reg. No.:

Awaiting Registration Call-Up

Incompatibility:

Anionic detergents, Natural Soaps, Strong oxidisers

INSTRUCTIONS FOR USE

Spray: Shake well before use. Spray **SteriTech Clinical Cleaner** onto the surface from a distance of about 20cm away - ensure that the surface is evenly covered by the foam. If not, spray more, or wipe the product over the uncovered area. Allow the minimum contact time of your required efficacy, as per the Product Efficacy table (below right) before rinsing, or using the surface again. Dispose as per the instructions provided on the label.

Wipes: Wipe the surface until visibly clean, allowing the minimum contact time of your required efficacy, as per the Product Testing table (below right) before using the surface again. Dispose as per the instructions provided on the label.



EFFICACY TESTING:



	Tested To:	Time: (mins)	Dilution:	Conditions:
Bactericidal Efficacy	SANS 51276 (EN 1276): Enterococcus hirae, Escherichia coli, Pseudomonas aeruginosa & Staphylococcus aureus SANS 53727 (EN 13727): Enterococcus hirae, Pseudomonas aeruginosa & Staphylococcus aureus	5	NEAT	Clean
Fungicidal & Yeasticidal Efficacy	SANS 51650 (EN 1650): Aspergillus brasiliensis & Candida albicans SANS 53624 (EN 13624): Aspergillus brasiliensis & Candida albicans	15	NEAT	Clean
Sporocidal Efficacy	SANS 53704 (EN 13704): Bacillus subtilis	15	NEAT	Clean
Mycobactericidal Efficacy	EN 14348 (SANS 54348): Mycobacterium terrae	60	NEAT	Clean
Corrosivity	SANS 636 (In-house): Non-corrosive up to 7 hours. Do not exceed a 7 hour exposure time.			
Dermal Irritation	OECD 439: Category 2 Irritant.			
Cytotoxicity	ISO 10993-5: Cytotoxic potential. Do not apply to broken skin, open wounds or cell-culture containers.			