



**SteriTech**  
SPECIALIST INFECTION PROTECTION

# REDUCE

## *SURFACE & EQUIPMENT DISINFECTANT WIPES*

### **ACTIVE INGREDIENTS**

0.072% Fifth generation QAC Blend;  
0.5% Preservative.

### **PRODUCT CHARACTERISTICS**

Vacuum-impregnated; Unfragranced;  
Non-corrosive; Non-toxic; Non-staining  
Non-caustic; Convenient

### **RECOMMENDED USES**

Surface disinfection.  
Equipment disinfection.  
Disinfection of high-touch items, especially  
where the use of a spray disinfectant is  
inadvisable.

### **PPE REQUIREMENTS**

**None strictly necessary**, however OSH  
requires the use of eye protection & good  
ventilation when mixing, handling or  
discarding chemical disinfectants.



For additional information about  
this product, or to give product  
feedback to the manufacturer,  
please scan this QR code.

## INSTRUCTIONS FOR USE

**To begin:** Remove the lid & remove the sealed inner bag. Ensure that the leading wipe is at the top, cut open the bag & place it back in the tub. Open the flap in the lid & pull the leading wipe slowly through the hole. Replace the lid onto the tub and the wipes are ready to use.

**When needed:** Simply pull out a wipe, tear it off on the perforation and use. Best results are achieved when used on a visibly clean, dry surface.

Allow the contact time of your required efficacy, (as per the efficacy table below right), before using the surface again. Repeat as often as required by your Infection Control protocol (suggested protocols available on request).

Always ensure that the flap is tightly shut after use to prevent drying out of the wipes.

**Wipe material is not efficiently biodegradable.**

## PRODUCT DESCRIPTION

**Available sizes:**

1000s Bucket (Wipe Size: 170mm x 170mm)  
 300s Bucket (Wipe Size: 170mm x 140mm)  
 100s Tub; 100s Refill (Wipe Size: 200mm x 140mm)

**Material Specification:**

20-30 gsm polypropylene / viscose blend

**Physical Form:**

White, non-woven wipe

**Odour:**

Mild chemical fragrance

**pH:**

6.5-7.5

**NRCS Registration:**

NRCS/8054/242613/248

**Incompatibility:**

Anionic detergents, Natural Soaps, Oxidisers



## EFFICACY TESTING:



	Tested To:	Time: (mins)	Dilution:	Conditions:
Bactericidal Efficacy	SANS 51276 (EN 1276): Enterococcus hirae, Escherichia coli, Pseudomonas aeruginosa & Staphylococcus aureus	5	NEAT	Clean
Fungicidal & Yeastical Efficacy	SANS 51650 (EN 1650): Aspergillus brasiliensis & Candida albicans	15	NEAT	Clean
Virucidal Efficacy	EN 13610 (SANS 53610): Lactococcus lactis P001 & Lactococcus lactis P008	15	NEAT	Clean
Sporocidal Efficacy	SANS 53704 (EN 13704): Bacillus subtilis	5	NEAT	Clean
Corrosivity	SANS 636 (In-house): Non-corrosive.			
Dermal Irritation	OECD 439: Non-irritant.			
Cytotoxicity	ISO 10993-5: Cytotoxic potential. Do not apply to broken skin, open wounds or cell-culture containers.			