



SteriTech
SPECIALIST INFECTION PROTECTION

SURFACE DISINFECTANT



GENERAL PURPOSE SURFACE DISINFECTANT

ACTIVE INGREDIENTS

0.1% Fifth generation QAC blend;
0.05% Biguanide; Surfactant

PRODUCT CHARACTERISTICS

Ready-to-Use; Fragranced;
Non-corrosive; Non-toxic; Non-staining
Non-caustic; Biodegradable

RECOMMENDED USES

Healthcare surface disinfection.
Medical Device disinfection.
NICU Incubator disinfection.
Eradication of mould, mildew and fungal
growth within healthcare facilities.



For additional information about
this product, or to give product
feedback to the manufacturer,
please scan this QR code.

PRODUCT DESCRIPTION

Available sizes:

25L, 5L, 1L Trigger Spray, 500mL Mist Spray, 100mL Mist Spray

Dilution Factor: Ready to use

Physical Form: Translucent green liquid

Odour: Mild lemon fragrance

pH: 9-9.5

Risk Classification: B / IIa

GMDN Code: 47631

GMDN Descriptor:

Medical device disinfection agent.

NRCS Registration:

NRCS/8054/242613/261

SAHPRA Medical Device Registration:

Awaiting Class Registration Call-up

Incompatibilities:

Anionic detergents, Natural Soaps, Strong oxidisers

PRODUCT PREVIOUSLY KNOWN AS STERISAFE AND STERI201

INSTRUCTIONS FOR USE

Thoroughly clean & rinse the surfaces prior to using the product. Then spray **SteriTech Surface Disinfectant** uniformly onto the surface, from a distance of about 30cm away. Allow the contact time of your required efficacy, as per the Efficacy Table below, or allow the surface to dry completely, naturally, before using the surface again. Discard expired or unused product as per the instructions on the label.

PPE REQUIREMENTS

None strictly necessary, however OSH requires the use of eye protection & good ventilation when mixing, handling or discarding chemical disinfectants.



EFFICACY TESTING:



	Tested To:	Time: (mins)	Dilution:	Conditions:
Bactericidal Efficacy	SANS 51276 (EN 1276): Enterococcus hirae, Escherichia coli, Pseudomonas aeruginosa & Staphylococcus aureus SANS 53727 (EN 13727): Enterococcus hirae, Pseudomonas aeruginosa & Staphylococcus aureus	5	NEAT	Clean
Fungicidal & Yeastidal Efficacy	SANS 51650 (EN 1560): Aspergillus brasiliensis & Candida albicans SANS 53624 (EN 13624): Aspergillus brasiliensis & Candida albicans	15	NEAT	Clean
Sporocidal Efficacy	SANS 53704 (EN 13704): Bacillus subtilis	15	NEAT	Clean
Virucidal Efficacy	EN 13610 (SANS 53610): Lactococcus lactis P001 & Lactococcus lactis P008 EN 14476 (SANS 54476): Human coronavirus	15	NEAT	Clean
Corrosivity	SANS 636 (In-house): Non-corrosive up to 7 hours. Do not exceed a 7 hour exposure time.			
Dermal Irritation	OECD 439: Non-irritant.			
Cytotoxicity	ISO 10993-5: Cytotoxic potential. Do not apply to broken skin, open wounds or cell-culture containers.			