



SteriTech
SPECIALIST INFECTION PROTECTION

DISINFECTING WIPES

*HEALTHCARE SURFACE &
EQUIPMENT DISINFECTANT
WIPES*

ACTIVE INGREDIENTS

0.072% Fifth generation QAC Blend;
0.5% Preservative.

PRODUCT CHARACTERISTICS

Vacuum-impregnated; Unfragranced;
Non-corrosive; Non-toxic; Non-staining
Non-caustic; Convenient

RECOMMENDED USES

Healthcare surface disinfection.
Medical Device disinfection.
Disinfection of high-touch items, especially
where the use of a spray disinfectant is
inadvisable.

PPE REQUIREMENTS

None strictly necessary, however OSH
requires the use of eye protection & good
ventilation when mixing, handling or
discarding chemical disinfectants.



PRODUCT DESCRIPTION

INSTRUCTIONS FOR USE

To begin: Remove the lid & remove the sealed inner bag. Ensure that the leading wipe is at the top, cut open the bag & place it back in the tub. Open the flap in the lid & pull the leading wipe slowly through the hole. Replace the lid onto the tub and the wipes are ready to use.

When needed: Simply pull out a wipe, tear it off on the perforation and use. Best results are achieved when used on a visibly clean, dry surface. Allow the contact time of your required efficacy, (as per the efficacy table below right), before using the surface again. Repeat as often as required by your Infection Control protocol (suggested protocols available on request).

Always ensure that the flap is tightly shut after use to prevent drying out of the wipes.

Wipe material is not efficiently biodegradable.

Available sizes:

1000s Bucket (Wipe Size: 170mm x 170mm)
300s Bucket (Wipe Size: 170mm x 140mm)
100s Tub; 100s Refill (Wipe Size: 200mm x 140mm)

Material Specification:

20-30 gsm polypropylene / viscose blend

Physical Form: White, non-woven wipe

Odour: Mild chemical fragrance

pH: 6.5-7.5

Risk Classification: B / IIa

GMDN Code: 58077

GMDN Descriptor:

Medical device cleaning/disinfection wipe

NRCS Registration:

NRCS/8054/242613/247

SAHPRA Medical Device Registration:

Awaiting Class Registration Call-up.

Incompatibility:

Anionic detergents, Natural Soaps, Oxidisers



EFFICACY TESTING:



	Tested To:	Time: (mins)	Dilution:	Conditions:
Bactericidal Efficacy	SANS 51276 (EN 1276): Enterococcus hirae, Escherichia coli, Pseudomonas aeruginosa & Staphylococcus aureus SANS 53727 (EN 13727): Enterococcus hirae, Pseudomonas aeruginosa & Staphylococcus aureus	5	NEAT	Clean
Fungicidal & Yeastical Efficacy	SANS 51650 (EN 1650): Aspergillus brasiliensis & Candida albicans SANS 53624 (EN 13624): Aspergillus brasiliensis & Candida albicans	15	NEAT	Clean
Virucidal Efficacy	EN 13610 (SANS 53610): Lactococcus lactis P001 & Lactococcus lactis P008 EN 14476 (SANS 54776): Human adenovirus	15 60	NEAT	Clean
Sporocidal Efficacy	SANS 53704 (EN 13704): Bacillus subtilis	5	NEAT	Clean
Corrosivity	SANS 636 (In-house): Non-corrosive.			
Dermal Irritation	OECD 439: Non-irritant.			
Cytotoxicity	ISO 10993-5: Cytotoxic potential. Do not apply to broken skin, open wounds or cell-culture containers.			