

ACCESSORIES



SMALL STERILISATION TRAY

500mL-1L sterilising tray for small instruments, complete with drip tray.

Dimensions (LxBxH):
280mm x 128mm x 90mm



MEDIUM STERILISATION TRAY

2L-4L sterilising tray for small-medium instruments, complete with drip tray.

Dimensions (LxBxH):
320mm x 160mm x 140mm



LARGE STERILISATION TRAY

4L-8L sterilising tray for large instruments & rigid scopes, complete with drip tray.

Dimensions (LxBxH):
640mm x 160mm x 140mm



KIDNEY STERILISATION TRAY

10L sterilising tray for flexibe endoscopes.

Approximate Dimensions (LxBxH):
900mm x 330mm x 200mm



M.E.C. TEST STRIPS (QAC)

Minimum Effective Concentration Test Strips for testing the active ingredient level of QAC-based products.



M.E.C. TEST STRIPS (GLUT)

Minimum Effective Concentration Test Strips for testing the active ingredient level of gluteraldehyde-based products.



M.E.C. TEST STRIPS (OPA)

Minimum Effective Concentration Test Strips for testing the active ingredient level of ortho-phthalaldehyde-based products.



WIRE FRAME WALL MOUNTED DISPENSER

For dispensing elbow pump (hospital) hand hygiene products in healthcare institutions. Bottle not included.



PUSH BUTTON WALL MOUNTED DISPENSER

Refillable dispenser for the manual dispensing of hand hygiene products. 1L reservoir capacity. Batteries not needed.



MANUAL WALL MOUNTED SACHET DISPENSER

Manual dispenser for hand hygiene products in sachet format. Batteries not needed. Reservoir not included.