



Chlorine Disinfectant Wipes >1000ppm

- High level disinfection
- Ready to use
- Rapid biocidal action
- Stabilised formulation

Medipal Chlorine Disinfectant Wipes are made with stabilised Sodium Hypochlorite NaClO. The wipes are less toxic and less corrosive than other chlorine releasing products whilst remaining effective against a broad spectrum of micro-organisms.

PASS EN16615



Medipal Chlorine Disinfectant Wipes are intended for the disinfection of non porous hard surfaces, equipment and non invasive medical devices within healthcare environments.

Our chlorine wipes have been independently tested to a range of European testing standards (EN) that are applicable to healthcare environments and are proven to be effective against viruses, bacteria and spores in minutes.

Broad spectrum efficacy

>1000ppm chlorine exhibits broad spectrum anti-microbial activity. Suitable for disinfecting isolation rooms, commodes and patient areas during outbreaks of infections.

Passes EN16615: 2015

Passes EN16615 in 1 minute. EN16615 is the first European standard aimed at testing the wipes efficacy rather than the extracted solution in the wipe.

Safer to use*

No need to dilute or activate the wipe
> 1000ppm chlorine dispensed onto the surface
No goggles or respiratory PPE required.

Fast acting

Full Virucidal*	1 minute
Bactericidal*	2 minute
<i>C. difficile</i>	2 minute
TB	1 minute

Please refer to the separate Efficacy Testing Summary Data sheet for the full results.

*As defined according to EN14885:2015

Widely compatible and safer to use*

Stabilised formulation of Sodium Hypochlorite making them safer to use.*

A lower pH of 9 compared with other chlorine products at pH 12.

Widely compatible, less corrosive and less damaging compared with other chlorine releasing products.

*Compared with other chlorine releasing products.



Packaging Type

Flow wrap

Code

S570110MPCE

Description

Medipal Chlorine Disinfectant Wipes 50

Count