



**PrimeSafe**  
Primary Infection Control

# ALCOHOL HAND RUB

(LIQUID & WIPES FORMAT)



## W.H.O. FORMULATION HAND SANITISER

### ACTIVE INGREDIENTS

**PrimeSafe Alcohol Hand Rub** contains  
80% v/v Ethanol; 1.45% Emollients;  
0.125% Hydrogen Peroxide;  
18.425% Inert Ingredients

### PRODUCT CHARACTERISTICS

WHO recommended formulation for alcohol based hand rubs. Non-staining, with a non-lingering alcohol odour.

### RECOMMENDED USES

Hand and skin disinfection.

As recommended by the Department of Health,  
U.S. C.D.C & W.H.O.

### PPE REQUIREMENTS

**None strictly necessary**, however OSH requires the use of eye protection & good ventilation when mixing, handling or discarding chemical disinfectants.



For additional information about this product, or to give product feedback to the manufacturer, please scan this QR code.



**THIS PRODUCT  
IS FLAMMABLE**

# INSTRUCTIONS FOR USE

## General Directions For Use (liquid):

Apply approximately 3mL (one pump) onto visibly clean hands and rub over all surfaces of the hands and wrists until dry. Do not rinse off.

## General Directions For Use (wipes):

To begin, remove the lid & ensure that the leading wipe is at the top of the tub. Open the flap in the lid & pull the leading wipe slowly through the hole. Replace the lid onto the tub and the wipes are ready to use. When needed, simply pull out a wipe, tear it off on the perforation and rub over all surfaces of visibly clean hands and wrists until dry. Do not rinse off. Allow to dry naturally before touching anything.

## Surgical Hand Preparation Directions (liquid only):

Apply a minimum of 3 applications (3 x 3mL) over a period of 3 minutes, to thoroughly pre-cleaned hands, ensuring that all areas of the hands and wrists are covered, and rub until dry. Then put on surgical gloves as per the procedure of the institution.

# PRODUCT DESCRIPTION

## Available sizes:

25L, 5L, 800mL Sachet, 1L Trigger Spray, 500mL Elbow Pump, 500mL Lotion pump, 500mL Trigger Spray, 500mL Mist Spray, 100mL Spray, 50mL Belt Clip Spray, 1000s Wipes Bucket (Wipe Size: 170mm x 170mm), 20s Wipe Purse Pack (Wipe Size: 170mm x 140mm), Individual Wipe Sachet (Wipe Size: 170mm x 140mm).

**Dilution Factor (Rub):** Ready to use

**Material Specification (Wipes):**

20-30 gsm polypropylene / viscose blend

**Physical Form (Rub):** Translucent liquid

**Physical Form (Wipe):** White, nonwoven wipe

**Odour:** Distinctive alcohol odour

**pH:** 5.0-7.5

**Risk Classification:** A / I

**GMDN Code:** 41550

**GMDN Descriptor:** Hand Disinfectant

**Incompatibility:**

Do not mix with other chemicals or products.



# EFFICACY TESTING:



	Tested To:	Time: (mins)	Dilution:	Conditions:
Bactericidal Efficacy	SANS 51276 (EN 1276): Enterococcus hirae, Escherichia coli, Pseudomonas aeruginosa & Staphylococcus aureus	1		
	SANS 53727 (EN 13727): Enterococcus hirae, Pseudomonas aeruginosa & Staphylococcus aureus	1	NEAT	Clean
	SANS 490:2020 (Alcohol-based Hand Rubs)	1		
	EN 1500 (Hygienic Hand Rubs)	3		
	EN 12791 (Surgical Hand Disinfection)	3		
Fungicidal & Yeastidal Efficacy	SANS 53624 (EN 13624): Aspergillus brasiliensis & Candida albicans	3	NEAT	Clean
Sporocidal Efficacy	SANS 53704 (EN 13704): Bacillus subtilis	1	NEAT	Clean
Virucidal Efficacy	EN 14476 (SANS 54476): Human adenovirus, Murine norovirus & Human coronavirus (EO)	3	NEAT	Clean
Irritation & Skin Corrosivity	SANS 490:2020 (Alcohol based Hand Rubs)	1		
	OECD 439: Episkin Model	3		
	EN 1500 (Hygienic Hand Rubs)	3	NEAT	Clean
	EN 12791 (Surgical Hand Disinfection)	3		