



PrimeSafe
Primary Infection Control

PURIFIED WATER



INSTRUCTIONS FOR USE

Use as required by the use instructions supplied with the product, chemical, equipment or medical device which necessitates the use of such water.

PPE REQUIREMENTS

None



SCAN THIS QR CODE FOR
ADDITIONAL INFORMATION OR
TO GIVE PRODUCT FEEDBACK TO
THE MANUFACTURER.

INTENDED USE

PrimeSafe Purified Water is water that has been filtered through a Reverse Osmosis process at 5-25 micro Siemens. It is intended for use in environments where water low in particulate matter and total dissolved solids is desirable.

INGREDIENTS

[100%] Reverse Osmosis filtered Water

PACKAGING VARIANTS

25L, 5L

PRODUCT CHARACTERISTICS

Recommended Dilution Factor: None

Physical Form: Clear liquid

Odour: None

pH 6.0-7.0

Incompatibilities: None

Environmental Concerns: None

Optional Accessories: None

REGULATORY INFORMATION

None applicable

ADDITIONAL INFORMATION

Product is not sterile.

