



HEXIRUB



HEALTHCARE SANITISER

ACTIVE INGREDIENTS

0.5% Chlorhexidine gluconate;
75% Ethanol (v/v), Emolient.

PRODUCT CHARACTERISTICS

Liquid format; Unfragranced;
Non-corrosive; Non-toxic; Non-staining
Non-caustic; Biodegradeable

RECOMMENDED USES

Healthcare skin disinfection.

PPE REQUIREMENTS

None strictly necessary, however OSH requires the use of eye protection & good ventilation when mixing, handling or discarding chemical disinfectants.



For additional information about this product, or to give product feedback to the manufacturer, please scan this QR code.



THIS PRODUCT IS FLAMMABLE

INSTRUCTIONS FOR USE

General Directions For Use:

Apply approximately 3mL of PrimeSafe HexiRub onto visibly clean hands and rub over all surfaces of the hands and wrists until dry. Do not rinse off. Or using a wipe, wipe all surfaces of the wrists and hands for at least 1 minute.

Surgical Hand Preparation Directions:

Apply a minimum of 2 applications (2 x 3mL) of PrimeSafe HexiRub, over a period of 3 minutes, to thoroughly pre-cleaned hands, ensuring that all areas of the hands and wrists are covered, and rub until dry. Then put on surgical gloves as per the procedure of the institution.

PRODUCT DESCRIPTION

Available sizes:

25L, 5L, 800mL Sachet, 500mL Elbow Pump, 500mL Lotion Pump, 500mL Mist Spray, 100mL Spray, 50mL Spray.

Dilution Factor:

Ready to use

Physical Form:

Clear liquid (previously pale blue)

Odour:

Distinctive alcohol odour

pH:

6.5-7.5

Risk Classification: B / IIa

GMDN Code: 41550

GMDN Descriptor: Hand Disinfectant

Incompatibility:

Do not mix with other chemicals



EFFICACY TESTING:



	Tested To:	Time: (mins)	Dilution:	Conditions:
Bactericidal Efficacy	SANS 51276 (EN 1276): Enterococcus hirae, Escherichia coli, Pseudomonas aeruginosa & Staphylococcus aureus	1		
	SANS 53727 (EN 13727): Enterococcus hirae, Pseudomonas aeruginosa & Staphylococcus aureus	3	NEAT	Clean
	SANS 490:2020 - Alcohol-based Hand Rubs	1		
	EN 1500 (SANS 5500) - Hygienic Hand Rubs EN 12791 (SANS 52791) - Surgical Hand Disinfection	3 3		
Fungicidal & Yeastocidal Efficacy	SANS 53624 (EN 13624): Aspergillus brasiliensis & Candida albicans	3	NEAT	Clean
Sporocidal Efficacy	SANS 53704 (EN 13704): Bacillus subtilis	1	NEAT	Clean
Dermal Irritation & Corrosivity	SANS 490:2020 - Alcohol-based Hand Rubs	1		
	OECD 439: Episkin model	3	NEAT	Clean
	EN 1500 (Hygienic Hand Rubs)	3		