



SteriTech
SPECIALIST INFECTION PROTECTION

DISINFECTING HANDWASH



ANTIMICROBIAL HAND CLEANSER

ACTIVE INGREDIENTS

1% Fifth generation QAC blend;
0.5% Biguanide; 8% Surfactant

PRODUCT CHARACTERISTICS

Ready-to-Use; Fragranced; Amphoteric;
Non-oily; Contains emollients;
Non-corrosive; Non-toxic; Non-staining
Non-caustic; Biodegradable

RECOMMENDED USES

Healthcare skin cleaning & disinfection.

PPE REQUIREMENTS

None strictly necessary, however OSH requires the use of eye protection & good ventilation when mixing, handling or discarding chemical disinfectants.

PRODUCT DESCRIPTION

Available sizes:

25L, 5L, 800mL Sachet,
500mL Elbow Pump; 500mL Lotion Pump

Dilution Factor:

Ready to use

Physical Form:

Translucent blue liquid

Odour:

Mild sweet cherry fragrance

pH: 6.5 - 7.5

Product Risk Classification: B / IIa

GMDN Code: 41551

GMDN Descriptor: Liquid Skin Soap

Incompatibilities:

Aldehydes and Acids

INSTRUCTIONS FOR USE

Thoroughly **wet hands first**, in order to achieve maximum foaming capacity. Apply one pump (approximately 3mL) of undiluted **SteriTech Disinfecting Handwash** onto wet hands. Work into a lather & spread over hands and wrists, ensuring a 1 minute contact time and using the appropriate hand washing technique.

For surgical hand disinfection allow a 3 minute contact time. Rinse well under running water. Dry hands thoroughly on a clean, disposable towel.



EFFICACY TESTING:



| | Tested To: | Time: (mins) | Dilution: | Conditions: |
|---------------------------------|---|--------------|-----------|-------------|
| Bactericidal Efficacy | SANS 51276 (EN 1276): Enterococcus hirae, Escherichia coli, Pseudomonas aeruginosa & Staphylococcus aureus | 1 | NEAT | Dirty |
| | SANS 53727 (EN 13727): Enterococcus hirae, Pseudomonas aeruginosa & Staphylococcus aureus | 1 | | |
| | EN 1499 (Hygienic Hand wash) | 1 | | |
| Fungicidal & Yeastidal Efficacy | SANS 51650 (EN 1650): Aspergillus brasiliensis & Candida albicans SANS 53624 (EN 13624): Aspergillus brasiliensis & Candida albicans | 3 | NEAT | Dirty |
| Sporocidal Efficacy | SANS 53704 (EN 13704): Bacillus subtilis | 3 | NEAT | Dirty |
| Dermal Irritation & Corrosivity | OECD 439 EN 1499 (Hygienic Hand wash) | 3 | NEAT | Dirty |